


## PRIMA FORM


## PRIMA PLOT


## KALLI HR KALLI VISCO



## ISKIO BM ISKIO VISCO

## KALLI DOS HR KALLI DOS VISCO

	PRIMA FORM	PRIMA PLOT	KALLI HR	KALLI VISCO	ISKIO BM	ISKIO VISCO	KALLI DOS HR	KALLI DOS VISCO
	20 < 90 kg 44 < 198 lb	20 < 100 kg 44 < 220 lb	-	20 < 100 kg 44 < 220 lb	30 < 150 kg 66 < 330 lb	30 < 150 kg 66 < 330 lb	-	-
	20 < 100 kg 44 < 220 lb	20 < 110 kg 44 < 242 lb	-	20 < 110 kg 44 < 242 lb	30 < 160 kg 66 < 352 lb	30 < 160 kg 66 < 352 lb	-	-
	37,5 - 50 cm	40 - 50 cm	38 - 50 cm	34 - 50 cm	38 - 50 cm	40 - 50 cm	43 / 48 cm	43 / 48 cm

\* Selon la référence / according to reference / nach der Referenz / según la referencia / volgens referentie / secondo il riferimento / i henhold til reference / de acordo com a referência / zgodnie z odniesieniem / viittauksen mukaan / عجز مدل ب حسب / prema referenci / Podľa referencie

	PRIMA FORM	PRIMA PLOT	KALLI HR	KALLI VISCO	ISKIO BM (Bi-matière)		ISKIO VISCO
(lxP)	(H = 7 cm)	(H = 7 cm)	(H = 8 cm)	(H = 8 cm)	(H = 8,5 cm)	(H = 10,5 cm)	(H = 10,5 cm)
34 x 36	-	-	-	17.17209	-	-	-
34 x 38	-	-	-	17.17212	-	-	-
36 x 36	-	-	-	17.17214	-	-	-
36 x 38	-	-	-	17.17216	-	-	-
37,5 x 38	17.17264	-	-	17.17234	-	-	-
38 x 38	-	-	04.234	-	-	-	-
38 x 40	-	-	-	-	17.17007	17.17107	-
38 x 42	-	-	-	17.17217	-	-	-
38 x 43	-	-	-	-	17.17112	17.17113	-
40 x 38	17.17265	17.17265P	-	-	-	-	-
40 x 40	17.17260	17.17260P	04.294	17.17294	17.17106	17.17108	17.17008
40 x 42	-	-	-	17.17218	-	-	-
40 x 43	-	-	-	-	17.17111	17.17114	-
40 x 45	-	-	-	17.17219	-	-	-
42 x 38	-	-	-	17.17228	-	-	-
42 x 40	-	-	-	17.17226	-	-	-
42 x 42	17.17266	-	04.295	17.17295/17.17296	-	-	-
42 x 45	-	-	-	17.17227	-	-	-
43 x 40	17.17270	17.17270P	-	-	-	-	-
43 x 43	-	-	-	-	17.17116	17.17118	17.17018
43 x 44	-	-	04.297	17.17267/17.17297	-	-	-
43 x 45	-	-	-	-	17.17032	17.17132	-
45 x 40	-	-	-	17.17225	-	-	-
45,5 x 40	17.17277	-	-	-	-	-	-
45,5 x 43	17.17275	17.17275P	-	-	-	-	-
45,5 x 45	17.17276	-	-	-	-	-	-
46 x 43	-	-	-	-	17.17126	17.17128	17.17028
46 x 45	-	-	04.235	17.17235	17.17034	17.17134	-
48 x 40	17.17279	-	-	17.17229	-	-	-
48 x 43	-	-	04.236	17.17236	-	-	-
48 x 45	-	-	-	17.17271	17.17133	17.17135	-
50 x 45	17.17280	17.17280P	04.299	17.17299	17.17136	17.17138	17.17038

	KALLI DOS HR		KALLI DOS VISCO	
(lxP)	(H = 44 cm)	(H = 46 cm)	(H = 44 cm)	(H = 46 cm)
43 x (5/10)	04.110	-	17.19110	-
48 x (5/10)	-	04.210	-	17.19210

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## 1. INDICATIONS

### Device Intended Use

- Prevention of pressure ulcers.
- Compensation for disability in the seated position (cushions and back).

### Indications

- Chair support for the prevention of pressure ulcers in patients who are at moderate to high-risk of developing pressure ulcers (as per clinical judgment and scales), with no pressure asymmetry and who are not-at risk for "forward-sliding" (KALLI VISCO, ISKIO VISCO, ISKIO BM, PRIMA FORM)
- Chair support for the prevention of pressure ulcers in patients who are at moderate to high-risk of developing pressure ulcers (as per clinical judgment and scales), with no pressure asymmetry, and who have a "forward-sliding" type of postural deficit (psychomotor regression, stroke, hypertonia of the hamstrings, hypertonia of the extensors of the spine) (PRIMA PLOT)
- Wheelchair seating support for patients with no risk of pressure ulcers (KALLI HR, back)
- Localised excessive pressure due to changes to the curvature of the spine. Lateral inclination of the trunk due to a mild postural deficit (back).

### Contra-indications

- Patients with frontal or sagittal stability difficulties.
- Established sacrococcygeal and/or ischial pressure ulcers.
- Patients whose weight is higher than that which the cushion can support.
- Risk of "forward sliding" (except for PRIMA PLOT).

### Target patient and user group

- Seated patients with reduced mobility with no pressure asymmetry and with no risk of sliding forward, who have a "moderate to high" risk for pressure ulcers (KALLI VISCO, ISKIO VISCO, ISKIO BM, PRIMA FORM).
- Seated patients with reduced mobility with no pressure asymmetry, presenting with a forward-sliding type postural deficit in the seated position, and who have a "moderate to high" risk for pressure ulcers (PRIMA PLOT).
- Seated patients with reduced mobility, with no risk of pressure ulcers (KALLI HR, back)

### Report undesirable side-effects

All serious incidents which are related to the device must be notified to the manufacturer and the competent authority of the member state in which the user and/or patient resides. Inform the competent authority if you consider or have reason to believe that the device presents a serious risk or that it is a fake product.

## 2. COMPOSITION OF THE MEDICAL DEVICE

### PRIMA PLOT

- 1 Viscoelastic polyurethane foam (80 kg/m<sup>3</sup>).
- 2 High resiliency abduction pad (40 kg/m<sup>3</sup>).
- 3 PHARMATEX cover, with anti-slide underside and transport handle.
- 4 Operating instructions (paper or electronic).

### KALLI VISCO, KALLI HR & PRIMA FORM

- 1 Viscoelastic polyurethane foam (80 kg/m<sup>3</sup>). (80 kg/m<sup>3</sup>) for KALLI VISCO and PRIMA FORM, or highly resilient mousse (37 kg/m<sup>3</sup>) for KALLI HR.
- 2 PHARMATEX cover, with anti-slide underside and transport handle. A JERSEY cover is also available for the KALLI HR.
- 3 Operating instructions (paper or electronic).

### ISKIO BM & ISKIO VISCO

- 1 Viscoelastic polyurethane foam (80 kg/m<sup>3</sup>).
- 2 Highly resistant foam insert (37 kg/m<sup>3</sup>) for ISKIO BM, or viscoelastic foam insert, with greater load-bearing capacity for ISKIO VISCO.
- 3 MAILLE 3D complete cover; or PHARMATEX cover, with anti-slide underside and transport handle.
- 4 Operating instructions (paper or electronic).

### BACK

- 1 Viscoelastic polyurethane foam (80 kg/m<sup>3</sup>). (80 kg/m<sup>3</sup>) for KALLI DOS VISCO, or highly resilient mousse (37 kg/m<sup>3</sup>) for KALLI DOS HR.
- 2 PHARMATEX cover (KALLI DOS VISCO). JERSEY cover and edges in PHARMATEX (KALLI DOS HR).
- 3 Operating instructions (paper or electronic).

## 3. CLINICAL BENEFIT, PERFORMANCE, MECHANISM OF ACTION

### Performance characteristics of the device

- Cushion for the prevention of pressure ulcers: The viscoelastic foam material allows for immersion of the pelvis, enveloping the areas in which bones stick out (ischium, coccyx, trochanter).
- Non-specific comfort cushions and back: Improved comfort in the seated position, that allows basic daily living tasks (bodily hygiene, transfers etc) to be performed.

### Expected clinical benefits

Maintains tissue oxygenation at the cutaneous and subcutaneous tissues that are in contact with the support, and/or contributes to the activities required for life support.

### Information for health professionals

Observe the condition of the patient's skin which is in contact with the cushion several times a day. Do not sit the patient down if they have a pressure ulcer at the ischial or coccygeal regions. Verify that the seat of the chair is in good condition. For wheelchairs, verify the height of the arm rests, the footrests etc. as well as its performance in terms of propulsion.

## 4. USAGE INSTRUCTIONS

### User training

User training must be carried out by persons who have been trained and approved by the relevant economic operators, especially in terms of safety and reporting of non-conformities.

### Device installation

The cushion is provided with the cover on (and with a spare cover in some cases) The product is ready to be installed. The anti-slide underside of the cushion must be in contact with the wheelchair, and the handle must be facing the front of the cushion. Sit the patient down.

### Preventive maintenance

- Perform an annual visual check of the condition of the foam.
- Perform an annual check of the external surface of the cover, by using a light source from inside in order to check for holes and/or tears.

### Cleaning and disinfection

- Between each patient.
- Wash the cover in a machine at a temperature of up to 90°C (PHARMATEX), or up to 40°C (MAILLE 3D).
- Do not iron the cover.
- Can be used with a bleach solution of up to 5,000 ppm (PHARMATEX). Cannot be used with a bleach solution (MAILLE 3D).

### Information regarding suitable procedures for re-use

- Each time that it is re-used, the product must be physically and bacteriologically cleaned by performing the procedures given in the technical manual (which can be downloaded at [www.winnicare.com](http://www.winnicare.com)).
- Replace the cover if the appearance of its surface changes.

## 5. WARNINGS AND PRECAUTIONS

### Required measures

- A cushion must be placed on a wheelchair in good condition.
- Ensure that the dimensions of the cushion are compatible with the seat of the wheelchair. The equipment must be used with its original protection.
- Check the position of the pelvis, the vertebral spine, the head, the armrests and the length of the leg rests and footrests after the patient has been seated on the cushion.
- A cushion by itself is not enough to prevent pressure ulcers/pressure wounds; other preventative measures are also required:
  - Change positions at least once every 2 hours.
  - Maintain skin hygiene and avoid maceration; do not put any other surface (sheet, towel etc.) between any part of the chair (back and seat) which would cause surface tension to appear, in order to maintain the efficacy of the cushion.
  - Once a day, observe, or ask someone else to observe, the condition of the skin.
  - In case of incontinence, regularly change the protective cover,
  - Regularly drink enough fluids.
- To ensure that the support is effective, it is important that limit excessive thicknesses between the body and the support (sheets, adult nappies etc). Mainly wear clothes that have a lot of cotton, if possible, without seams in the support region.
- Ensure that no foreign bodies are present, such as tubing, crumbs, fatty substances etc.

➤ Important: please talk with your prescriber in the month following the purchase of the cushion.

➤ Enhanced monitoring for patients who are at risk of falls (patient profile: disabled, agitated, not lucid).

### Circumstances under which the user must talk with a healthcare professional

- If one of these measures cannot be implemented, you must inform your doctor or nurse as soon as possible.
- Inform your doctor or nurse as soon as possible of any abnormal events such as fever, pain, reddening or whitening of the weight-bearing areas (ischium, coccyx, trochanter, hip, shoulder blade).

### Precautions for use

- Bone and/or muscular and/or skin trauma in the pelvic and thigh or spinal regions, that is not stabilised and/or painful, and that is in direct contact with the support device.
- Instability in the seated position.
- The first few days following surgery for a pressure ulcer/wound (skin graft or flap).
- Please also verify the condition of the exposed part of the skin each time that care is provided.

### Warnings

(Re)evaluate the patient's risk of entrapment in the non-moving parts of the wheelchair equipped with the therapeutic cushion.

### Lifetime information

The device may be reused within the initial performance maintenance limits for the duration of its life cycle. Its estimated lifetime is of 4 years.

### Accessories

Only protective accessories for the cushion supplied by the manufacturer can ensure the safety and performance of the complete product.

## 6. STORAGE/ HANDLING/DISPOSAL

### Specific conditions for storage and/or handling

Store away from light, heat, and humidity.

### Disposal of the product

Do not throw the product out into non-dedicated locations in nature. Follow the recycling channels that are in place in your country.

Pressures (max) in mmHg before ageing, using the maximum validated weight

	Max. user weight	Max. user pressure Buttocks
PRIMA FORM PRIMA PLOT	100 kg	126
KALLI VISCO	130 kg	169
ISKIO BM	160 kg	104
ISKIO VISCO	150 kg	104

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Siotec AG, Sägestrasse 75, CH-3098 KÖlnig
Mangar International Limited (WINNCARE group)
Cubro Limited. 149 Taurikura Drive, Tauriko,