



# Standing slings

**ergo**  
slings

## User manual

- Standing sling, 1-belt
- Standing sling, 2-belts

## Important information

Read this user manual before using the product.

This user manual is aimed at care personnel, family members and therapists in connection with the use and maintenance of slings made by ErgoSling.





# Contents

1. PURPOSE .....	4
2. SAFETY INSTRUCTIONS .....	4
3. CHECK BEFORE USE.....	5
4. LIFTING.....	6
5. STANDING SLINGS.....	7
6. MAINTENANCE .....	8
6.1 Washing instructions .....	8
6.2 Storage.....	8
7. INSPECTION .....	8
8. QUALITY AND THE ENVIRONMENT.....	9
8.1 Disposal.....	9
8.2 Quality .....	9
8.3 Symbols and abbreviations.....	9

## 1. PURPOSE

This range of slings from **ErgoSlings** are all standing sling that are intended for use for supporting people with reduced mobility in a transfer situation.

Transfer can be achieved in combination with a stand-aid lifter.

- These slings can be used at hospitals, in the context of home care, at institutions and in private homes.
- The slings can be used for transfers from a wheelchair/chair or from a chair to a bathing and toilet chair etc.
- All models of the sling are available in various sizes. See the specifications for the various sling models in this user manual.

**NOTE:** All sling models described in this manual. Illustrations and descriptions are applicable to all versions, unless a separate description is provided.

**NOTE:** **ErgoSlings** must not be used for any purposes other than those described in this user manual.

## 2. SAFETY INSTRUCTIONS



### Important – Warning!

***Read the user manual and these safety instructions thoroughly before using the sling:***

- Do not exceed the sling's maximum lifting capacity.
- Use the correct size and model of sling for the user.
- Make sure the sling has no visible damage or signs of wear. In particular check that the seams are intact. The hoist straps must not be frayed, and they must not be knotted.
- For hygiene reasons, the sling must be clean. See the washing instructions in section 5.
- Qualified personnel must always be present when attaching the sling and during lifting.
- It must be ensured that hazardous situations do not arise when using any other equipment, *e.g.*: the sling straps becoming caught in the wheels of a wheelchair, or similar.
- It must be ensured that neither the helper nor the user can be injured while using the sling.
- The user must never be left alone in the sling.
- The helper is responsible for ensuring that the sling is used properly and is also responsible for ensuring that safety is maintained. See section - 3. *Check before use.*

### 3. CHECK BEFORE USE

The helper is responsible for ensuring that the sling is checked before use, and that instructions for use have been provided.

**Check the following points before use:**

- 1) Check that the sling used is the correct model and is the right size, so that it is right for the lifting situation and the user.

**Note:** All the slings have a coloured loop that indicates the size and the centre line of the sling.



**Size chart** (all measurements are approx.)

	XS	S	M	L	XL	XXL
	62 – 80 cm	66 – 90 cm	80 – 100 cm	100 – 126 cm	110 – 132 cm	120 – 140 cm
Colour						

- 2) Check that the fabric of the sling and the straps are intact. There must not be any wear or damaged seams.
- 3) Check that the sling is clean.
- 4) If in any doubt about how to use the sling, Winnicare Nordic can be contacted for advice and guidance.

## 4. LIFTING

**Note:** The lifting situation must always be planned, and it must be ensured that the user feels secure when being lifted in the sling.

- 1) Place the sling behind the user's back with the colored loop by the label indicating the centerline. The centerline must be aligned with the user's spine.
- 2) The non-slip piece must face the user's back.
- 3) The standing sling with 2-buckles has a velcro-piece that must be tightly closed around the user.



- 4) Click the buckle / buckles together and tighten the belt.
- 5) Attach the straps on the lifting arms hooks. See instructions on page 7.
- 6) Now lift the user carefully, when the straps are stretched the sail can be straightened so that it sits correctly.



## Attach the sling to the spreader bar



WinnCare Nordic excludes any liability if the sling is not correctly attached to the hook on the spreader bar.

It is always the responsibility of the helper to ensure that the correct sling is used and that it is suitable for the need/purpose

- 1) Attach the sling strap to the open hook.
- 2) Pulling the strap down causes the lock to close with a “click”.



**Note:** The lifting hook is designed in such a way that the user's weight automatically keeps the hook locked.

To unhook the strap, pulling the strap up causes the hook to open.



Care should be taken to avoid the user's fingers becoming caught in the hook.

## 5. STANDING SLINGS

### Standing sling with one buckle



### Standing sling with two buckles



### Use:

- Standing sling is for users who can stand by themselves during the transfer, but needs assistance getting up from sitting position or assistance to sit down.
- Standing sling are used for users with standing function and a good upper body control.

## Specification:

Standing sling with one buckle	
SWL	250 kg
Material	Polyester
Padding	Yes
Sling is available in sizes	S - XXL
Standing sling with two buckles	
SWL	250 kg
Material	Polyester
Padding	Yes
Sling is available in sizes	XS - XL

## 6. MAINTENANCE

### 6.1 Washing instructions

**ErgoSlings** can be machine-washed. We recommend that the buckles on the slings are closed and washed in a laundry bag to protect the buckles.

**Note:** Do not use softener.

Standing sling can be tumble-dried at a low temperature.

Follow the washing instructions as indicated by symbols and printed instructions on the sling label.

A key to the symbols can be found under section 7.3 - *Symbols and abbreviations*

### 6.2 Storage

If the sling is not being used, it should be packed in its original packaging and stored in a clean and dry place at a temperature between 0-40°C.

## 7. INSPECTION

**ErgoSlings** should be inspected at least once every six months in accordance with EN/ISO 10535, Annex B, which recommends that inspection should be carried out by a qualified person who is familiar with use of the sling.

The legal owner of the product is responsible for ensuring that inspection is carried out. An inspection schedule can be requested from Winncare Nordic

### Warranty

There is a two-year guarantee on materials and manufacture defects, on condition that there is compliance with safety and washing instructions.



## 8. QUALITY AND THE ENVIRONMENT

**ErgoSlings** have an expected lifetime of 5 - 10 years assuming normal daily use. This also assumes that the recommended inspection and maintenance is carried out.

### 8.1 Disposal

**ErgoSlings** may be disposed of together with flammable materials at the recycling centre, or they can be returned to Winncare Nordic.

### 8.2 Quality

Winncare Nordic is certified in accordance with ISO 9001 and ISO 13485 quality standards.

These certificates mean that the product and the company meet applicable international standards for quality management and standards for traceability for medical products.

**ErgoSlings** are classified as Medical equipment class 1 product.

### 8.3 Symbols and abbreviations



Machine-wash at the indicated temperature.



Do not iron.



Do not use bleach.



Do not dry-clean.



Tumble-dry at low temperature.



See the user manual.



The product can be disposed of with general waste.

**SWL**

Safe working load = (user + accessories)



The product is CE-marked.





---

**Fabricator:**

WinnCare Nordic Ltd  
Taarnborgvej 12 C  
4220 Korsoer  
Denmark

Phone: + 45 70 27 37 20  
Fax: + 45 70 27 37 19  
Email: [salg@winnCare.dk](mailto:salg@winnCare.dk)  
[www.winnCare.dk](http://www.winnCare.dk)

93/42/EØF CLASS I

DS/EN 12182



910-00769UK