

Luna Spreader bars




2-point spreader bar
For bayonet or strap
Available in 360/480/580 mm



Straight 2-point spreader bar
Available in 360/480/580 mm



Room-to-room 2-point spreader bar
For Luna 275 kg
Available in 360/480/580 mm



3-point spreader bar
For bayonet or strap



4-point spreader bar
For bayonet or strap



3-point powered spreader bar
For Luna 275 kg



Spreader bars for Luna Twin, 400 kg
Løftebøjle og fleksible side-løftebøjler
Available in 360/480/580 mm



2-point spreader bar on fixed Luna

Key standard functions and highlights


- Extensive range for specific and general use
- Patented sling hooks on all spreader bars
- Easily interchangeable to allow for multiple use
- For both standard and bariatric lifting situations
- Modular to achieve multi flexibility
- Patented safety system on hooks
- Extensive sling range available

Item Numbers

2-point spreader bar	360 mm	480 mm	580 mm
For bayonet	015-02129	015-02126	015-02128
For strap	015-02108	015-02110	015-01968
2-point straight spreader bar	360 mm	480 mm	580 mm
	015-02113	015-02114	015-02116
Room-to-room spreader bar	360 mm	480 mm	580 mm
	015-02117	015-02118	015-02119
3-point spreader bar			
For bayonet	015-02106		
For strap	015-02104		
4-point spreader bar			
For bayonet	015-02102		
For strap	015-02100		
3-point powered spreader bar			
For Luna 275 kg	015-02112		
Luna Twin (400 kg) spreader bars	360 mm	480 mm	580 mm
2-point spreader bar	015-02206	015-02204	015-02202
Flexible site-spreader bar	015-02036	015-02037	015-02038



Spreader bar for bayonet on Luna.
Spring loaded locking mechanism.
To mount and dismount: Push in. Turn. Release.



Sling attachment: a minimum of downward pressure locks the sling hook; pull upwards to unlock sling hook.

Please consult the technical manual for further specifications.
WinnCare Nordic reserves the right to changes without notice.
Product photos may include accessories. Please consult your dealer for further information.

CE MDD93/42/EEC
Class I

date 04.03.2021 - ver. 3

WinnCare Nordic ApS
Tårnborgevej 12C
4220 Korsør
Danmark
tlf. +45 70 27 37 20
e-mail: info@winnCare.dk

www.winnCare.dk

WINCARE
N o r d i c