



## TOILET SLING

Perfect half sling, allowing good contact between the user and the carer

### USAGE

- For transfers to and from the toilet.
- For users with a good control of the head and torso.
- For transfers from a sitting or semi-sitting position.
- Used for transfers in mobile and stationary hoists.

### FEATURES

- The sling has a large opening around lower body, which makes toileting and changing of clothes easier.
- The sling provides good support under the thigh and the back of pelvis to the shoulder blade.
- The sling ensures a good upright posture that facilitates the positioning of the user.
- Toilet sling has split padded leg straps.
- Location pocket for better sling fitting when sitting in a wheelchair.
- Straps on both sides of the sling to help the carer moving the user.



### SPECIFICATIONS

SWL: 275 kg

Material: Polyester

Sizes: XS-S-M-L-XL-XXL

### ITEM NUMBERS

Toilet sling - polyester

	Item number
<b>XS</b>	009-09543
<b>S</b>	009-09544
<b>M</b>	009-09545
<b>L</b>	009-09546
<b>XL</b>	009-09547
<b>XXL</b>	009-09548

## SLING SIZE GUIDE

The sling has a coloured strap on upper back with size and midline mark, which ensures correct placement of the sling.

Colour code is presented on the size chart.

The image below shows a red colour code, which means that the sling is in size S (Small):



	XS	S	M	L	XL	XXL
<b>Y</b> (hip measurement)	37 cm	42 cm	48 cm	52 cm	68 cm	78 cm
<b>X</b> (measurement from neck to lower back)	50 cm	55 cm	66 cm	75 cm	75 cm	75 cm



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. WinnCare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.

## FOR MORE INFORMATION VISIT

[winnCare.com](http://winnCare.com)

E-mail: [info@winnCare.dk](mailto:info@winnCare.dk)

Phone: +45 70273720

WIN  CARE