



UNIVERSAL BASIC with head support

Supports thighs and body. Great seating comfort.

USAGE

- For users with lack of stability of torso. The sling with head support is used for patients, who need extra support around the head.
- For transfers from a sitting or lying position.
- For transfers in mobile and stationary hoists.

FEATURES

- The sling is transferring the user from a wheelchair/chair or from lying position from a bed or bath.
- The sling provides perfect support around the thighs and the high back, supporting the upper body and giving the user a feeling of comfort and safety.
- The sling ensures a good upright posture that facilitates the positioning of the user.
- The mesh sling is suitable for bathing.
- Location pocket for better sling fitting when sitting in a wheelchair.
- Straps on both sides of the sling to help the carer moving the user.



SPECIFICATIONS

SWL: 275 kg

Material:

Polyester and polyester mesh.

Sizes: XS-S-M-L-XL-XXL

ITEM NUMBERS:

Universal Basic with head support, polyester:

	Item no.
XS	009-09823
S	009-09824
M	009-09825
L	009-09826
XL	009-09827
XXL	009-09828

Universal Basic with head support, polyester mesh:


	Item no.
XS	009-09833
S	009-09834
M	009-09835
L	009-09836
XL	009-09837
XXL	009-09838

SLING SIZE GUIDE

The sling has a coloured strap on upper back with size and midline mark, which ensures correct placement of the sling.

Colour code is presented on the size chart.

The image below shows a yellow colour code, which means that the sling is in size M (Medium):



	XS	S	M	L	XL	XXL
Y (hip measurement)	37 cm	42 cm	48 cm	52 cm	68 cm	78 cm
Z (measurement from head to lower back)	75 cm	80 cm	91 cm	100 cm	100 cm	100 cm



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. Winncare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.



FOR MORE INFORMATION VISIT

winncare.com

E-mail: info@winncare.dk

Phone: +45 70273720

WIN  CARE